

Protector® Special Application Laboratory Fume Hoods





Protecting you.

Unique applications require unique fume hoods. Our special application fume hoods are here to protect you. From radiochemicals to harsh acids, we have you covered.

Mining precious metals. Etching electrical wafers for everyday electronics. Creating components for rocket fuel. From radiochemicals to harsh acids, there's a Labconco hood tailored to fit your requirements.

Not sure what you need? Our ventilation specialists will analyze your application and recommend the right hood for your laboratory.

We've been taking care of customers since 1925. Let us put our experience to work for you.

Table of Contents

Protector® Stainless Steel Radioisotope Laboratory Hoods		Specifications	6
Features & Benefits	3	Ordering Information & Dimensional Data	7
Specifications	4	Protector® PVC Perchloric Acid & Acid Digestion Laboratory Hoods	
Ordering Information & Dimensional Data	5	Features & Benefits	8
Protector® Stainless Steel Perchloric Acid Laboratory Hoods		Specifications	9
Features & Benefits	3	Ordering Information & Dimensional Data	10
		Washdown System Specifications	11

Protector® Stainless Steel Radioisotope & Perchloric Acid Laboratory Hoods

Features, Benefits & Selection Guide

By-pass airflow design
for stable face velocities.
Restricted by-pass increases
air velocity above the sash.

**Stainless steel baffle and
liner with integral work
surface** for a seamless work
area. Perchloric Acid Hoods
also have a drainage trough.

Built-in washdown system*
for removing hazardous
perchlorates from behind the
baffle. It consists of
pre-plumbed internal piping,
spray nozzles and upper
left side-mounted remote
control knob.

**Optional color-coded
service fixtures** for gas, air,
water, vacuum and other
services have remote
controls for use regardless of
the sash position. Left corner
post is factory-prepared for
up to three service fixtures,
right for up to four (eight
fixtures total per hood).
**Additional service fixtures
require holes drilled in the
liner at the factory or on site.**

❖ **Large unobstructed
sightline** of 37.5" high from
the work surface to the
header panel allows taller
users comfortable viewing
while standing.

❖ **Clean-Sweep™ sash handle and tracks**

Openings on the sash handle bleed air into the hood chamber and away from the operator's breathing zone. Slots on the sash tracks enhance airflow.

Bright interior lighting

High-efficiency, instant
start, T8 LED lights are
located outside the hood
interior for corrosion-
resistance and easy
replacement.

❖ **Pre-wired electrical components**

LED lights and switches are
factory-wired to the hood's
single point junction box.
Each hood is factory-
prepared for up to four
electrical duplexes and
an airflow monitor.

Cord-Keeper™ slots on
both sides of the air foil
allow the sash to close
completely when electrical
cords from equipment inside
the hood are plugged in.
Cords are kept out of the
way of the operator.



❖ **Eco-Foil™ air foil with Clean-Sweep™ openings** reduces energy consumption by 7-10% compared to flat air foils. Aerodynamic curve allows air to sweep the work surface for maximum containment. Air foil is durable powder-coated stainless steel.

	Liner Material	Washdown System	Applications	Pages
Radioisotope Hood	Type 304 stainless steel withstands high temperatures	—	Radiochemicals Nuclear medicine Micro electronics Contact Labconco for use with lead bricks or shields	4-5
Perchloric Acid Hood	Type 316 stainless steel withstands high temperatures and is less chemically reactive than Type 304 stainless steel	Includes valve, drain trough and nozzles (see page 11)	Perchloric acid	6-7

*Washdown system and drainage trough on Perchloric Acid Hoods only.

❖ Labconco exclusive feature

Protector® Stainless Steel Radioisotope Laboratory Hoods

Specifications



All models feature:

- By-pass airflow design and by-pass block
- Glacier white powder-coated steel exterior
- ✱ Eco-Foil air foil with aerodynamic Clean-Sweep airflow openings*
- ✱ Cord-Keeper slots on left and right side of air foil
- Type 304 stainless steel liner with radiused corners, integral work surface and pre-set baffle(s)
- ✱ Powder-coated aluminum sash handle with Clean-Sweep openings and Clean-Sweep slots on the powder-coated stainless steel sash tracks*
- Tempered safety glass vertical-rising sash with cable and pulley
- Removable front and side panels, and front service panels for access to plumbing and electrical wiring
- Pre-wired T8 LED lighting with vapor-proof design and ADA-compliant light and blower switches
- Black enamel-coated hardboard supporting surface
- Type 304 stainless steel, 11.0" OD exhaust connection(s)

Standards conformance & regulations:

- CFR 29, Part 1910
- SEFA 1
- NFPA 45
- ASTM E84-09C
- ASHRAE 110
- ANSI Z9.5
- UL 61010-1
- CAN/CSA C22.2 No. 61010.1
- UL 1805

✱ Labconco exclusive feature

- CE Conformity Marking (230V models)
- SEFA 8

Fixtured models may feature:

- Two pre-plumbed service fixtures with forged brass valves, lower right side with brass tubing for gas and lower left side with copper tubing for air. Components for converting either or both fixtures to cold water and vacuum are provided. **Inlet tubing is not provided**
- One pre-wired GFCI electrical duplex receptacle on lower right side and, on 8' models only, one additional pre-wired GFCI electrical duplex receptacle on lower left side

Required accessories (not included):

- Dedicated remote blower
- Ductwork. Type 304 stainless steel duct is recommended. Contact Labconco for source recommendations
- Base cabinet or stand

Optional accessories for on-site installation include:

- Service Fixture Kits
- Guardian Airflow Monitor Kits
- Electrical Duplex Receptacle Kits
- Sash Stop Kits
- Ceiling Enclosure and Rear Finish Panel Kits



Total Exhaust CFM and Static Pressure

Face Velocity (fpm)	Airflow Volumetric Rate (CFM) @ Static Pressure (inches of water)							
	4' Hood		5' Hood		6' Hood		8' Hood	
	CFM	s.p.	CFM	s.p.	CFM	s.p.	CFM	s.p.
Sash @ Full Open (28")								
125	905	0.39"	1195	0.52"	1475	0.78"	2050	0.44"
100	725	0.25"	955	0.33"	1180	0.50"	1640	0.28"
Sash @ 62.5% Open (18")								
125	565	0.15"	745	0.20"	920	0.31"	1280	0.17"
100	455	0.10"	595	0.13"	740	0.20"	1025	0.11"

*U.S. Patent No. 6,461,233

♿ Heights of switches, electrical receptacles and service fixtures meet requirements of Americans with Disabilities Act (ADA).

Protector® Stainless Steel Radioisotope Laboratory Hoods

Ordering Information & Dimensional Data

Catalog Number Configurator:

Use this key to configure the nine digit catalog number to order your Protector Radioisotope Laboratory Hood.

1 2 0 0 0

Select the **WIDTH** of your fume hood. This is the fourth digit of your catalog number. Shipping weight is also noted for 31.7" models. Add 65 lbs. (29 kg) for 37.7" models. Add 10 lbs. (5 kg) for fixtured models.


4 = 4' (122 cm)	6 = 6' (183 cm)
720 lbs. (327 kg)	835 lbs. (379 kg)
5 = 5' (152 cm)	8 = 8' (244 cm)
775 lbs. (351 kg)	955 lbs. (433 kg)

Select the **DEPTH** of your fume hood. This is the fifth digit of your catalog number.


0 = 31.7" (80.5 cm)
1 = 37.7" (95.8 cm)

Select the **ELECTRICAL REQUIREMENTS, SERVICE FIXTURES** and **GFCI ELECTRICAL DUPLEX RECEPTACLE** combination you desire. These are the eighth and ninth digits of your catalog number.


Electrical Requirements	No Service Fixtures	Two Service Fixtures	Two Service Fixtures & GFCI Duplex*
100-115V, 50/60 Hz, 10A	00	01	02
208-230V, 50/60 Hz, 5A	20	21	—



TOP



FRONT



SIDE

Visit labconco.com for detailed AutoCAD drawings.

*Hoods with GFCI electrical duplex are rated at 20 amps. 8' hoods have two GFCI electrical duplex receptacles, one mounted on each side, rated at 20 amps each.

**8' model has two exhaust collars, 24.0" on each side of centerline.

Protector® Stainless Steel Perchloric Acid Laboratory Hoods

Features & Benefits



All models feature:

- By-pass airflow design and by-pass block
- Glacier white powder-coated steel exterior
- ✱ Eco-Foil air foil with aerodynamic Clean-Sweep airflow openings*
- ✱ Cord-Keeper slots on left and right side of air foil
- Type 316 stainless steel liner with radiused corners, integral work surface, drainage trough and pre-set baffle(s)
- ✱ Powder-coated aluminum sash handle with Clean-Sweep openings and Clean-Sweep slots on the powder-coated stainless steel sash tracks*
- Tempered safety glass vertical-rising sash with cable and pulley
- Removable front and side panels, and front service panels for access to plumbing and electrical wiring
- Pre-wired T8 LED lighting with vapor-proof design and ADA-compliant light and blower switches
- Built-in washdown system with internal piping and spray nozzles (see specifications on page 11)
- Black enamel-coated hardboard supporting surface
- Type 316 stainless steel, 11.0" OD exhaust connection(s)

Standards conformance & regulations:

- CFR 29, Part 1910
- SEFA 1
- NFPA 45
- ASTM E84-09C
- ASHRAE 110
- ANSI Z9.5

- UL 61010-1
- CAN/CSA C22.2 No. 61010.1
- UL 1805
- CE Conformity Marking (230V models)
- SEFA 8

Fixtured models may feature:

- Two pre-plumbed service fixtures with forged brass valves, lower right side with brass tubing for gas and lower left side with copper tubing for air. Components for converting either or both fixtures to cold water and vacuum are provided. **Inlet tubing is not provided**
- One pre-wired GFCI electrical duplex receptacle on lower right side and, on 8' models, one additional pre-wired GFCI electrical duplex receptacle on lower left side

Required accessories (not included):

- Dedicated remote PVC blower
- Wash rings and other ductwork
- Base cabinet or stand

Optional accessories for on-site installation include:

- Service Fixture Kits
- Washdown Valve Kit
- Guardian Airflow Monitor Kits
- Electrical Duplex Receptacle Kits
- Sash Stop Kits
- Ceiling Enclosure and Rear Finish Panel Kits



Total Exhaust CFM and Static Pressure

Face Velocity (fpm)	Airflow Volumetric Rate (CFM) @ Static Pressure (inches of water)							
	4' Hood		5' Hood		6' Hood		8' Hood	
	CFM	s.p.	CFM	s.p.	CFM	s.p.	CFM	s.p.
Sash @ Full Open (28")								
125	905	0.39"	1195	0.52"	1475	0.78"	2050	0.44"
100	725	0.25"	955	0.33"	1180	0.50"	1640	0.28"
Sash @ 62.5% Open (18")								
125	565	0.15"	745	0.20"	920	0.31"	1280	0.17"
100	455	0.10"	595	0.13"	740	0.20"	1025	0.11"

*U.S. Patent No. 6,461,233

✱ Labconco exclusive feature

♿ Heights of switches, electrical receptacles and service fixtures meet requirements of Americans with Disabilities Act (ADA).

Protector® Stainless Steel Perchloric Acid Laboratory Hoods

Ordering Information & Dimensional Data

Catalog Number Configurator:

Use this key to configure the nine digit catalog number to order your Protector Stainless Steel Perchloric Acid Laboratory Hood.

1 3 0 0 0

Select the **WIDTH** of your fume hood. This is the fourth digit of your catalog number. Shipping weight is also noted for 31.7" models. Add 65 lbs. (29 kg) for 37.7" models. Add 10 lbs. (5 kg) for fixtured models.


4 = 4' (122 cm)	6 = 6' (183 cm)
770 lbs. (349 kg)	885 lbs. (401 kg)
5 = 5' (152 cm)	8 = 8' (244 cm)
825 lbs. (374 kg)	1000 lbs. (454 kg)

Select the **DEPTH** of your fume hood. This is the fifth digit of your catalog number.


0 = 31.7" (80.5 cm)
1 = 37.7" (95.8 cm)

Select the **ELECTRICAL REQUIREMENTS, SERVICE FIXTURES** and **GFCI ELECTRICAL DUPLEX RECEPTACLE** combination you desire. These are the eighth and ninth digits of your catalog number.

Electrical Requirements	No Service Fixtures	Two Service Fixtures	Two Service Fixtures & GFCI Duplex*
100-115V, 50/60 Hz, 10A	00	01	02
208-230V, 50/60 Hz, 5A	20	21	—




TOP



FRONT

	A	B	C
4' Hood	48.0" (121.9 cm)	38.25" (97.2 cm)	24.0" (61.0 cm)
5' Hood	60.0" (152.4 cm)	50.25" (97.2 cm)	30.0" (76.2 cm)
6' Hood	72.0" (182.9 cm)	62.25" (158.1 cm)	36.0" (91.4 cm)
8' Hood	96.0" (243.8 cm)	86.25" (219.1 cm)	24.0" (61.0 cm)**



SIDE

Visit labconco.com for detailed AutoCAD drawings.

*Hoods with GFCI electrical duplex are rated at 20 amps. 8' hoods have two GFCI electrical duplex receptacles, one mounted on each side, rated at 20 amps each.

**8' model has two exhaust collars, 24.0" on each side of centerline.

Protector® PVC Perchloric Acid & Acid Digestion Laboratory Hoods

Features, Benefits & Selection Guide

By-pass airflow design
for stable face velocities.
Restricted by-pass increases
air velocity above the sash.

**Corrosion-resistant baffle
and liner of Type 1
unplasticized PVC with
integral work surface** is
crevice-free and withstands
surface temperatures up to
140° F (60° C).

Built-in washdown system
facilitates the removal of
hazardous perchlorates from
behind the baffle. It consists
of pre-plumbed internal
piping, spray nozzles and
upper left side-mounted
remote control knob.

**Optional color-coded
service fixtures** for gas, air,
water, vacuum and other
services have remote
controls for use regardless of
the sash position. Left corner
post is factory-prepared for
up to three service fixtures,
right for up to four (eight
fixtures total per hood).
Additional service fixtures
require holes drilled in the
liner at the factory or on site.

❖ **Clean-Sweep™ sash handle and tracks**
Openings on the sash handle bleed air into the hood chamber and away from
the operator's breathing zone. Slots on the sash tracks enhance airflow.


Bright interior lighting
High-efficiency, instant start,
T8 LED lights are located
outside the hood interior for
corrosion-resistance and
easy replacement.

**Pre-wired electrical
components**
❖ LED lights and switches are
factory-wired to the hood's
single point junction box.
Each hood is factory-
prepared for up to four
electrical duplexes and an
airflow monitor.

**Large unobstructed
sightline** of 37.5" high from
the work surface to the
header panel allows taller
users comfortable viewing
while standing.

Cord-Keeper™ slots on
❖ both sides of the air foil
allow the sash to close
completely when electrical
cords from equipment inside
the hood are plugged in.
Cords are kept out of the
way of the operator.

❖ **Eco-Foil™ air foil with Clean-Sweep™ openings** reduces energy consumption
by 7-10% compared to flat air foils. Aerodynamic curve allows air to sweep
the work surface for maximum containment. Air foil is durable powder-coated
stainless steel.



	Liner Material	Sash Material	Applications	Pages
Perchloric Acid Hood	Type 1 PVC	Tempered safety glass	Perchloric acid and other inorganic chemicals that do not involve high temperatures	9-10
Acid Digestion Hood	Type 1 PVC	Lexan® (polycarbonate)	Hydrofluoric acid and other mineral acids that do not involve high temperatures Not for use with perchloric acid	9-10

*Lexan is a registered trademark of SABIC Innovative Plastics.

❖ Labconco exclusive feature

Protector® PVC Perchloric Acid & Acid Digestion Laboratory Hoods

Specifications



All models feature:

- By-pass airflow design and by-pass block
- Glacier white powder-coated steel exterior
- * Eco-Foil air foil with aerodynamic Clean-Sweep airflow openings*
- * Cord-Keeper slots on left and right side of air foil
- Heat-welded Type 1 unplasticized PVC liner with integral work surface, drainage trough and pre-set baffle(s). PVC withstands maximum continuous operating temperature of 140° F (60° C) and has a flame spread less than 25 per ASTM E84
- * Powder-coated aluminum sash handle with Clean-Sweep openings and Clean-Sweep slots on the powder-coated stainless steel sash tracks*
- Removable front and side panels, and front access panels for access to plumbing and electrical wiring
- Pre-wired T8 LED lighting with vapor-proof design and ADA-compliant light and blower switches
- Built-in washdown system with internal piping and spray nozzles (see page 11 for specifications)
- Black enamel-coated hardboard supporting surface
- PVC, 12.75" OD exhaust connection(s)

Standards conformance & regulations:

- CFR 29, Part 1910
- SEFA 1
- NFPA 45
- ASTM E84-09C
- ASHRAE 110
- ANSI Z9.5

*U.S. Patent No. 6,461,233

* Labconco exclusive feature

- UL 61010-1
- CAN/CSA C22.2 No. 61010.1
- UL 1805
- CE Conformity Marking (230V models)
- SEFA 8

⚠ **WARNING:** Cancer - P65Warnings.ca.gov (California only)

Protector PVC Perchloric Acid Hoods feature:

- 3/16" thick tempered safety glass vertical-rising sash with cable and pulley. Designed for perchloric acid use

Protector PVC Acid Digestion Hoods feature:

- 1/4" thick Lexan polycarbonate vertical-rising sash with cable and pulley. Designed for hydrofluoric acid use

Fixtured models may feature:

- Two pre-plumbed service fixtures with forged brass valves, lower right side with brass tubing and plastic serrated hose connector for gas and lower left side with copper tubing and gray PVC swivel gooseneck faucet for cold water. Components for converting the lower right fixture to air or vacuum are provided.
Inlet tubing is not provided
- One pre-wired GFCI electrical duplex receptacle on lower right side and, on 8' models, one additional pre-wired GFCI electrical duplex receptacle on lower left side

Required accessories (not included):

- Dedicated remote PVC blower
- Wash rings and other ductwork
- Base cabinet or stand

Optional accessories for on-site installation include:

- Service Fixture Kits
- Guardian Airflow Monitor Kits
- Electrical Duplex Receptacle Kits
- Sash Stop Kits
- Ceiling Enclosure and Rear Finish Panel Kits



Total Exhaust CFM and Static Pressure

Face Velocity (fpm)	Airflow Volumetric Rate (CFM) @ Static Pressure (inches of water)							
	4' Hood CFM	s.p.	5' Hood CFM	s.p.	6' Hood CFM	s.p.	8' Hood CFM	s.p.
Sash @ Full Open (28")								
125	905	0.27"	1195	0.44"	1475	0.59"	2050	0.35"
100	725	0.17"	955	0.28"	1180	0.38"	1640	0.22"
Sash @ 62.5% Open (18")								
125	565	0.11"	745	0.17"	920	0.23"	1280	0.13"
100	455	0.07"	595	0.11"	740	0.15"	1025	0.09"

⚠ Heights of switches, electrical receptacles and service fixtures meet requirements of Americans with Disabilities Act (ADA).

Protector® PVC Perchloric Acid & Acid Digestion Laboratory Hoods

Ordering Information & Dimensional Data

Catalog Number Configurator:

Use this key to configure the nine digit catalog number to order your Protector PVC Laboratory Hood.

1 4 1 0 0

Select the **APPLICATION TYPE** of your fume hood. This is the third digit of your catalog number.

0 = Perchloric Acid (with safety glass sash)

1 = Acid Digestion (with Lexan sash)

Select the **WIDTH** of your fume hood. This is the fourth digit of your catalog number. Shipping weight is also noted. Add 10 lbs. (5 kg) for fixtured models.

4 = 4' (122 cm)

735 lbs. (333 kg)

6 = 6' (183 cm)

850 lbs. (386 kg)

5 = 5' (152 cm)


790 lbs. (358 kg)

8 = 8' (244 cm)


970 lbs. (440 kg)

Select the **ELECTRICAL REQUIREMENTS, SERVICE FIXTURES** and **GFCI ELECTRICAL DUPLEX RECEPTACLE** combination you desire. These are the eighth and ninth digits of your catalog number.


Electrical Requirements	No Service Fixtures	Two Service Fixtures*	Two Service Fixtures & GFCI Duplex**
100-115V, 50/60 Hz, 10A	00	01	02
208-230V, 50/60 Hz, 5A	20	21	—



TOP



FRONT



SIDE

Visit labconco.com for detailed AutoCAD drawings.

*Left side service fixture is a Swivel Gooseneck Faucet. Right side service fixture is serrated hose connector.

**Hoods with GFCI electrical duplex are rated at 20 amps. 8' hoods have two GFCI electrical duplex receptacles, one mounted on each side, rated at 20 amps each.

†8' model has two exhaust collars, 24.0" on each side of centerline.


Protector® Stainless Steel & PVC Laboratory Hoods

Washdown System Specifications

Rinsing the fume hood prevents the accumulation of reactive perchlorate salts.

Hoods with washdown feature:

- Built-in washdown system with internal piping and spray nozzles with maximum flow rate from 0.5 - 1.0 GPM
- Include upper left side-mounted remote control fixture with forged brass valve with maximum flow rate of 3 GPM and maximum working pressure ranging from 10 to 40 psi depending on number of hood's spray nozzles
- Valve is capable of supplying water to hood's spray nozzles and one wash ring (wash rings not included) on 4', 5' and 6' models. Valve is capable of supplying water to hood's spray nozzles on 8' models
- Washdown fitting located on top of hood to facilitate connection to external wash rings (wash rings not included)



Washdown System Flow Rates for Protector Stainless Steel Laboratory Hoods

	Number of Spray Nozzles	Flow @ 10 psi	Flow @ 40 psi
4' Hood	3	0.6 gpm	1.2 gpm
5' Hood	4	0.8 gpm	1.6 gpm
6' Hood	4	0.8 gpm	1.6 gpm
8' Hood	6	1.2 gpm	2.4 gpm

Washdown System Flow Rates for Protector PVC Laboratory Hoods

	Number of Spray Nozzles	Flow @ 10 psi	Flow @ 40 psi
4' Hood	3	1.56 gpm	3.0 gpm
5' Hood	4	2.08 gpm	4.0 gpm*
6' Hood	4	2.08 gpm	4.0 gpm*
8' Hood	6	3.12 gpm	6.0 gpm*

*Each washdown valve can only produce up to 3.5 gallons per minute at 40 psi. These operational points require an additional Washdown Valve Kit (9808310).

Everything you need to succeed.

Ventilated Enclosures & Exhausters

Fume Hoods & Enclosures
Filtered Fume Hoods
Ductless (Carbon-Filtered)
Enclosures
Biological Safety Cabinets &
Enclosures
Clean Benches
PCR Enclosures
Balance, Bulk Powder &
Equipment Enclosures
Nanotechnology Enclosures
Laboratory Animal
Research Stations
Non-ventilated Enclosures
Remote Blowers & Exhausters

Sample Preparation Equipment

Freeze Dryers
Nitrogen Blow Down
Dry Evaporators
Vacuum Dry Evaporators
Vacuum Concentrators
Cold Traps
Vacuum Pumps

Glassware Washers

Water Purification Systems

Glove Boxes & Gas Purifiers

Controlled Atmosphere
Glove Boxes
Filtered Glove Boxes
Combination Glove Boxes

Forensic Enclosures

Kjeldahl, Fat & Crude Fiber Apparatus

Other Products

Laboratory Carts,
Tables & Benches
Blood Drawing Chairs
Vacuum Desiccator
Base Stands, Base Cabinets &
Work Surfaces
Filters & Prefilters
Vacuum Pumps
Freeze Dryer Accessories
Ventilation Product Accessories



**Visit labconco.com to request information
for all Labconco products.**

LABCONCO CORPORATION

8811 Prospect Avenue
Kansas City, MO 64132

(800) 821-5525 | (816) 333-8811

labconco.com

©2018 LABCONCO CORPORATION

Printed in the U.S.A.

Product design subject to change without notice.

2-57-05/24/24

LABCONCO®