



TOP VIEW



FRONT VIEW



SIDE VIEW

GENERAL DESCRIPTION

3' Paramount Ductless Enclosure, Tall-360 Visibility
Model: 6963302

- Exterior Superstructure
- Tempered Glass
- Dry powder epoxy-coated aluminum
- Interior Superstructure
- Dry powder epoxy-coated Aluminum and Steel Slides and Back
- 1/4" Thick Tempered Safety Glass
- Filters
- Accommodates 4 Filters
- Gas-assisted filter clamping mechanism

Electrical Fixtures

- Pre-wired Fluorescent Light with Timers

Control Panel / LCD Display

- Face Velocity
- Filter Status and Types
- Touchpad Control of Enclosure Operation
- Set Filter Life and Check Timers
- Programmable Alarm Set Points
- Event timer

Wiring

- Wired at Factory
- 8/2 Power Cord with Plug

SASH

- Gas-Assisted
- Pivots to full open for cleaning and material transfer
- 1/4" Tempered Safety Glass
- 12" Operating Height
- Epoxy Coated Aluminum Frame

Baffles

- 1/4" Tempered Safety Glass, Pivots Downward for Cleaning
- Air Foil
- Extruded aluminum, Clean Sweep Airfoil

Motor

- ECM Impeller, Smart-Flow System to Adjust Motor Speed

Sensor

- Safety First Vapor Sensor
- Primary and Secondary Back-Up Sensors

Utility Port

- (2) 1.5" diameter utility ports with plugs
- Assembly and Crating
- Pre-Assembled
- Politized

Exhaust

- Exhaust Shroud, recirculates filtered air back into room

FUME HOODS for:



Labconco Corporation
8811 Prospect Avenue
Kansas City, Mo 64132-2696
(800) 821-5525
www.labconco.com

No. _____ Date _____

Date: 4-7-2009

Drawing # 6963302

Drawn By: Luke Savage

Checked By:

Approved By:

Sheet No.